

A Comparison of IRMS and CRDS for the Detection of Economic Adulteration of Lemon Juice and Honey



In the last several years, economically motivated adulteration (EMA) of foods has received increased attention. Two popular targets are lemon juice and honey. For both of these types of adulteration analyses, the FCC is evaluating the use of Cavity Ring **Down Spectroscopy (CRDS) to establish its suitability** in comparison to the well established technique Isotope Ratio Mass Spectrometry (IRMS)

- Addition of C4 material will alter natural C3 ¹³C/¹²C ratio, indicative of adulteration.
- Iemon Juice adulterated with exogenous citric acid and/or low cost sweeteners \rightarrow C4 plant source
- Honey adulterated with low cost sweeteners \rightarrow C4 plant sugars.
- Lemon juice samples are imported, part of an FDA assignment.
- Honey samples bought from local markets.
- **Basic Procedure: Citric acid precipitate from lemon** juice or bulk honey and protein precipitate \rightarrow combusted in IRMS and CRDS for determination of ¹³C/¹²C (converted to per mil basis (δ^{13} C(‰)).
- \succ δ^{13} C citric acid less negative than -23‰ is indicative of lemon juice adulteration.
- $\succ \delta^{13}C_{\text{Protein}} \delta^{13}C_{\text{Honey}}$ less negative than -1‰ is indicative of honey adulteration.



- **PICARRO Combustion Module by Costech with** G2121i Isotopic CO₂ analyzer.
- \succ At start, N₂ carrier gas is switched to a volume of O₂ depending on size and composition of the sample.
- > Sample in tin capsule is dropped into combustion reactor, combusts at 1700-1800°C, breaks down into N₂, CO₂, H₂O and SO₂.
- Cu wires absorb excess O₂, SO₂ adsorbed by oxidation catalysts water adsorbed by Water Trap.
- \succ N2 + CO₂ \rightarrow GC column \rightarrow Liaison Interface module \rightarrow CO₂ analyzer $\rightarrow \delta^{13}$ C/¹²C measurement. CO₂ analyzers use time based, optical absorption spectroscopy



Samples are introduced into combustion oven at 980 °C, produces gases, then into reduction oven at 650 °C, get reduced. The gases pass through chemical trap and moisture trap, into the separating column, where CO₂ is separated from N₂ and other gases, and into the detector via electron impact ionization. Ions with m/z 44, 45 and 46 are collected in separate Faraday Cups and a complex algorithm calculates the ¹³C/¹²C ratio.

 δ <u>Sample</u> $= \left(\frac{R_{s}}{R_{s}}\right)$ Reference





- signature of the analyte gas.
- long.
- ring down time.

Madhavi Mantha, John R Urban and Kevin M. Kubachka US FDA, Forensic Chemistry Center, Cincinnati, OH ASITA June 2016, Philadelphia, PA

$$\frac{\text{Sam} - \text{R}_{\text{Ref}}}{\text{R}_{\text{Ref}}} \int 10^3 = \left(\frac{\text{R}_{\text{Sam}}}{\text{R}_{\text{Ref}}} - 1\right) 10^3$$

CRDS Schematic and Operation

Light from a tunable semiconductor diode laser is directed into an optical resonator cavity containing the analyte gas. When the optical frequency matches the resonance frequency of the cavity, energy builds up in the cavity. When the build-up is complete, the laser is shut off. The energy decays from the cavity exponentially in time, or "rings down," with a characteristic decay time. This energy decay is measured, as a function of time, on a photodiode. The ring down time is measured at several different wavelengths as the laser is tuned across the molecular

CRDS is a measurement of time not of absorbance. When the laser is at a wavelength where the gas in the cavity is strongly absorbing, the ring down time is short; when the wavelength is such that the gas does not absorb, the ring down time is

The concentration is proportional to the difference in these

Method Development and Validation for lemon juice

- Doner L., "Carbon Isotope Ratio in Natural and Synthetic Citric Acid as Indicators of Lemon Juice Adulteration, J. Agric. Food Chem., 22(1985) 770-772
- Method development to optimize calcium citrate precipitation Varied amounts of Ca(OH)₂ Heating time
- Changing the pH with NaOH and NH₄OH
- Measureable CO₂ background signal From carbonate contamination in Ca(OH), reagent present in calcium citrate precipitate.
- On achieving dissatisfactory results, focus was shifted to Using Ca(NO₃)₂ and CaCl₂ Varying the conditions as mentioned above
- Finally using 2 mL of 3M CaCl₂ solution at pH greater than 8.5 and heating at 60°C for 2 hours, produced calcium citrate precipitate, free of carbon contamination

Validation

Two commercially available citric acids Three locally purchased lemon juice from concentrates Freshly squeezed locally purchased lemons

Reference Materials (RM) and Quality Control

- > NIST RM 8542 (Sucrose, $\delta^{13}C_{VPDB} = -10.45 \pm 0.07\%$)
- > NIST RM 8573 (L-glutamic Acid, $\delta^{13}C_{VPDB}$ = -26.39±0.09‰)
- Acetanilide used to condition reactors and verify measurement performance for δ^{13} C.
- \succ RMs are used to normalize instrument δ^{13} C to sample based $\delta^{13}C_{VPDB}$, using Excel for calculations

Advantages and Disadvantages

EA-IRMS

- Well established technique with good linearity and precision
- High maintenance
- High degree of technical knowledge needed for operation
- Very expensive
- Used for ¹³C/¹²C, ²H/¹H, ¹⁸O/¹⁶O, ³⁴S/³²S and ¹⁵N/¹⁴N
- > CM-CRDS
- Relatively inexpensive
- Robust
- Simplified analysis technique
- Relatively new technique with no officially defined methods.
- Only works for ¹³C/¹²C analysis



	IRMS			CR	
Sample Number	Honey	Protein	Protein - Honey	Honey	Protein
	$\delta^{13}C[VPDB]$	$\delta^{13}C[VPDB]$	δ ¹³ C[VPDB]	$\delta^{13}C[VPDB]$	δ ¹³ C[VPD
15-1	-26.19	-26.00	0.19	-26.27	-26.33
16-1	-26.18	-26.19	-0.01	-26.39	-26.00
16-2	-25.45	-25.86	-0.40	-25.48	-25.80
16-3	-26.38	-26.16	0.22	-26.46	-26.70
16-4	-25.66	-25.87	-0.21	-25.90	-26.40
16-5	-25.19	-25.51	-0.33	-25.11	-25.56
Suspect-A Sub 1	-25.78	-26.92	-1.15	-25.86	-27.02
Suspect-A Sub 4	-25.71	-26.84	-1.13	-25.74	-26.89
Suspect-A Sub 8	-25.74	-26.95	-1.21	-25.75	-27.01
Suspect-B Sub 1	-25.64	-25.81	-0.17	-25.57	-26.00
Suspect-C Sub 1	-24.71	-25.57	-0.86	-24.74	-25.68

- Eleven honey samples were analyzed for Bulk Honey and Protein analysis and results were compared.
- **Eight of the eleven samples had** $(\delta^{13}C_{Protein} \delta^{13}C_{Honev})$ less negative than -1‰, indicative of no honey adulteration. Three samples were adulterated.
- \succ The difference in $\delta^{13}C[VPDB]$ (‰) values between IRMS and CRDS for Bulk Honey ranged from -0.07‰ to + 0.24‰.
- \succ The difference in $\delta^{13}C[VPDB]$ (‰) values between IRMS and CRDS for Bulk Honey ranged from -0.06‰ to + 0.54‰.





-1.26

-0.43

- \succ The average difference in $\delta^{13}C[VPDB]$ (‰) values between IRMS and CRDS for bulk honey was 0.28‰.
- The above results indicate that the new CM-CRDS is quite suitable for adulteration analysis when applied to lemon juice and honey samples.

Future Work

> To explore the CM-CRDS instrument on more honey samples and extend to other matrices like maple syrup and other fruit juices.

Acknowledgements