

Easier Analyzing of Stable Carbon Isotopes of Dissolved Inorganic Carbon in Water Samples or Solid Carbonate Samples: New Picarro Caddy Continuous Flow Interface + AutoMate Prep Device Pairing

PICARRO
OS23E-1469

Juan Carlos Guerrero¹ (jcguerrero@picarro.com), Jan Woźniak¹, Jason H. Curtis^{2,3}

(1) Picarro, Inc., Santa Clara, CA, USA, (2) University of Florida, Department of Geological Sciences, Gainesville, FL, USA, (3) AutoMate FX, Inc., Bushnell, FL, USA

1. BACKGROUND

Abstract

The Caddy Continuous Flow Interface has previously been used to connect available solid and liquid bulk sample preparation instruments—including the Picarro Combustion Module—to Picarro analyzers for high-precision carbon isotope ($\delta^{13}\text{C}$) measurements.


Picarro is happy to announce the AutoMate Prep Device from AutoMate FX as the Caddy's newest pairing. The Caddy's fully automated, low-cost, and simple operation will take advantage of the newest AutoMate Prep Device design improvements to bring together increased data quality and reduced interface complexity.

Here we present some results from tests done with the new Caddy + AutoMate Prep Device pairing, focused on assessing the system's precision at the lower end of our specification range for CO_2 concentrations. The results show the same or at times, improved precision levels and reduced memory effects otherwise associated with the previous interface device. All while making the pairing simpler to set up, easier to operate, and lower in cost.

2. OBJECTIVES

1. Determine how well the Caddy + AutoMate Prep Device pairing operates with CO_2 concentrations towards the lower end of our specification range (<380 ppm CO_2) by viewing precision values.
2. Assess memory effect between samples to confirm reduced memory effects seen during previous tests.
3. Confirm successful "proof-of-concept" tests previously conducted at Picarro HQ.

3. PICARRO CADDY + AUTOMATE PREP DEVICE PAIRING



The Picarro Caddy Continuous Flow Interface + AutoMate Prep Device system is composed of three components:


1. Picarro isotopic analyzer (e.g., G2201-i or G2131-i)
2. Picarro Caddy Continuous Flow Interface
3. AutoMate Prep Device manufactured by AutoMate FX, Inc.

4. LOW CO_2 CONCENTRATION STUDY

Overview

Dissolved inorganic carbon (DIC) solutions were made by mixing different weighted amounts of either sodium carbonate (Na_2CO_3) or sodium bicarbonate (NaHCO_3) to achieve CO_2 concentration values at or near the lower end of the analyzer's specification range. The AutoMate Prep Device adds phosphoric acid to samples to liberate CO_2 , which flows to the Picarro for $\delta^{13}\text{C}$ - CO_2 analysis. These ranged between 360 – 395 ppm CO_2 for Na_2CO_3 solutions and between 330 – 360 ppm CO_2 for NaHCO_3 solutions.

Results – Na_2CO_3 Samples




Each sample was analyzed three times. Below is a table showing the averaged results and precision per Na_2CO_3 solution sample.

| Sample | Avg $\delta^{13}\text{C}$ (‰) | Precision (‰, 1 σ) | Avg [CO_2] (ppm) |
|-----------|-------------------------------|----------------------------|-----------------------------|
| Sample #1 | -2.747 | 0.075 | 377.3 |
| Sample #2 | -2.833 | 0.107 | 385.1 |
| Sample #3 | -2.749 | 0.035 | 370.1 |
| Sample #4 | -2.834 | 0.096 | 387.9 |
| Sample #5 | -2.871 | 0.148 | 388.8 |
| Sample #6 | -2.892 | 0.070 | 384.3 |

Precision between 0.035 – 0.148‰

Results – NaHCO_3 Samples



Each sample was analyzed three times, except for Sample #1 due to an error (analyzed twice). Below is a table showing the averaged results and precision per NaHCO_3 solution sample.

| Sample | Avg $\delta^{13}\text{C}$ (‰) | Precision (‰, 1 σ) | Avg [CO_2] (ppm) |
|-----------|-------------------------------|----------------------------|-----------------------------|
| Sample #1 | 0.2685 | 0.0134 | 343.2 |
| Sample #2 | 0.3797 | 0.1045 | 347.3 |
| Sample #3 | 0.3320 | 0.0394 | 339.3 |
| Sample #4 | 0.3867 | 0.0180 | 342.6 |


Precision between 0.0134 – 0.1045‰

5. ASSESSING MEMORY EFFECT

Overview

Dissolved inorganic carbon (DIC) solutions were made by mixing weighted amounts of either sodium carbonate (Na_2CO_3) or sodium bicarbonate (NaHCO_3) to achieve CO_2 concentration values between 430 – 480 ppm CO_2 for Na_2CO_3 and NaHCO_3 solutions. The AutoMate Prep Device adds phosphoric acid to samples to liberate CO_2 , which flows to the Picarro for $\delta^{13}\text{C}$ - CO_2 analysis. Each sample was analyzed four times in an interleaved sample series.

Results



Below is a table showing the averaged results and precision per solution sample.

| Sample | Avg $\delta^{13}\text{C}$ (‰) | Precision (‰, 1 σ) | Avg [CO_2] (ppm) |
|---------------------------------------|-------------------------------|----------------------------|-----------------------------|
| Na_2CO_3 – 1.86 mg #1 | -2.331 | 0.086 | 449.5 |
| NaHCO_3 – 0.9 mg #1 | 0.681 | 0.134 | 439.0 |
| Na_2CO_3 – 1.86 mg #2 | -2.318 | 0.062 | 466.7 |
| NaHCO_3 – 0.9 mg #2 | 0.749 | 0.135 | 446.8 |
| Na_2CO_3 – 1.86 mg #3 | -2.274 | 0.101 | 457.9 |
| NaHCO_3 – 0.9 mg #3 | 0.719 | 0.051 | 453.3 |
| Na_2CO_3 – 1.86 mg #4 | -2.302 | 0.080 | 452.8 |

Precision between 0.062 – 0.101 for Na_2CO_3 samples

Precision between 0.051 – 0.135 for NaHCO_3 samples

6. CONCLUSIONS



- The Picarro Caddy + AutoMate Prep Device pairing had precision values always better than **0.2%**, and on average **~0.08%**.
- The Picarro Caddy helps **remove memory effect**, otherwise associated with the sample bags inside the previous interface device.
- The Picarro Caddy + AutoMate Prep Device pairing is an **overall simpler set up** compared to the previous interface device, while also being **easier to operate** and **lower in cost**.