

Improved throughput for $\delta^{18}\text{O}$ and δD measurements of water with Cavity Ring-Down Spectroscopy.



Improved throughput for $\delta^{18}\text{O}$ and δD measurements of water with Cavity Ring-Down Spectroscopy (CRDS)

Keren Drori^{1*}, Magdalena E. G. Hofmann¹, Lin, Zhiwei¹, Bent, Jonathan D.¹, and Jan Woźniak¹

¹Picarro, Inc., Santa Clara, USA

* Contact: kdrori@picarro.com

INTRODUCTION

- Oxygen ($^{18}\text{O}/^{16}\text{O}$) and deuterium (D/H) isotopes are a widespread tool to trace physical and chemical processes in hydrology and bio-geosciences. Precision and throughput are key parameters for water isotope analysis.
- Picarro presents two new methodologies for the L2130-i CRDS water isotope analyzer that allow the user to increase the throughput without compromising data quality.
- The Picarro **Express Method** distinguishes between a memory reduction stage and a sample analysis stage and allows the user to measure up to 50 samples per day while maintaining the excellent precision of CRDS (i.e., 0.01‰ for $\delta^{18}\text{O}$ and 0.05‰ for δD). This corresponds to doubling the throughput compared to the standard Picarro methodology.
- The Picarro **Survey Method** makes use of ultrafast injections and sorts the samples by their measured isotopic values, enabling a powerful new strategy to reduce memory effects.

INSTRUMENTATION

Picarro L2130-i water isotope analyzer

- Analyzer is coupled to an autosampler and vaporizer (see picture top right)
- Current specification: 15 permeg (SD of averages of 6 injections, i.e., 1 h measurement)
 - The current 15 permeg precision specification is limited by the reproducibility of the water injections, not by the precision of the analyzer.
 - An improved injection procedure (see box below) allows significantly higher throughput (see results).

METHODS

Methodological Improvements in Water Isotope Analysis




Fig. 1 Standard vs. Express mode. (b) Standard water injection procedure for a Picarro water isotope analyzer (L2140-i / L2130-i) coupled to an autosampler and vaporizer. All injections show the same peak shape. The first three injections are discarded due to memory effects. (a) The new injection procedure (Express) distinguishes between memory removal and sample injections. The six memory removal injections are done within a few minutes so that the water vapor peak reaches up to 70,000 ppm. The actual sample measurements are still done at 20,000 ppm H_2O but with a shorter integration time on the peak and a faster transfer from the vaporizer to the analyzer.

RESULTS

Sample Sorting Can Reduce Memory Effect




Fig. 2 Survey Method. Water vapor is sticky; the memory is hard to eliminate, even with multiple wet flushes. 18 different water samples were measured using Standard (○) or Survey (●) modes. The difference between two consecutive samples is significant in the original order for both ^{18}O (a) and δD (b). The best memory reduction efficiency is achieved when the samples are sorted by the δD values surveyed in the original sample order. Both the real $\delta^{18}\text{O}$ (c) and δD (d) values follow the new sample order nicely in a monotonic way, which indicates a good sorting.

Increase Throughput Without Compromise on Data Quality




Fig.3 Express Method. WICO 2020 samples were measured using the new survey (□), Express (x), and Standard (○) modes. Both ^{18}O (a) and δD (b) measurements showed excellent correlation between the Standard and the Express mode, suggesting similar precision. Measurements done using Survey mode correlate with the true value and can be used to re-sort samples.

CONCLUSIONS

The new water injection procedure allows increased measurement precision and sample throughput.

- Using Express mode can increase sample throughput by a factor of 2 and reach the same precision for $\delta^{18}\text{O}$ and δD as Standard mode.
- Using Survey mode allows for re-sorting of samples to reduce memory effect resulting in better precision.

REFERENCES

- [1] A. Pierchala, K. Rozanski, M. Dulinski, Z. Gorczyca, M. Marzec, R. Czub, High-precision measurements of $\delta^2\text{H}$, $\delta^{18}\text{O}$ and $\delta^{17}\text{O}$ in water with the aid of cavity ring-down laser spectroscopy. *Isotopes Environ. Health Stud.* **55**, 290–307 (2019).
- [2] E. Barkan, B. Luz, High precision measurements of $^{17}\text{O}/^{16}\text{O}$ and $^{18}\text{O}/^{16}\text{O}$ ratios in H_2O . *Rapid Commun. Mass Spectrom.* **19**, 3737–3742 (2005).

PICARRO