
PICARRO G2132-*i* δ¹³C High Precision Isotopic CH₄ CRDS Analyzer

High precision, real-time carbon isotope information from methane at ambient concentrations

- Isotopic measurements with superb precision and stability for critical methane source identification
- Direct measurement of CH₄ no sample drying or preparation required
- Field deployable: < 26 kg, meets shock and vibration tests
- Less hassle less calibration, less maintenance, no consumables
- Equipped with ChemDetect[™] to sense contaminants in the samples



The core of the G2132-*i* is Picarro's unique Cavity Ring-Down Spectroscopy (CRDS) technique, a time-based measurement that uses a laser to quantify spectral features of gas phase molecules in a small optical cavity, which has an effective laser path length of up to 20 kilometers. The G2132-*i* also includes a patented, high precision wavelength monitor, which maintains absolute spectral position, thereby ensuring accurate peak quantification of the absorption features for ¹²CH⁴ and ¹³CH⁴. Like all Picarros, the G2132-*i* maintains consistently high precision and stability due to its unique controls for temperature and pressure which would otherwise dramatically impact spectral measurements.



Excellent measurements of CH_4 are only possible if the influences of other molecules in the sample are quantified. The G2132-i makes precise measurements of CO_2 and H_2O concentrations to allow the cross-influence of these species to be quantified and corrected while providing more insight into the sample. The analyzer comes equipped with ChemDetectTM, a new layer of analysis that inspects the recorded spectra to find indications of contaminating species.

Portable, Rugged and Easy to Use: The G2132-*i* is compact, weighing less than 26 kg, which makes it easy to take to the field. It can be running within minutes out of the box and will continue to run for months without user interaction. The G2132-*i* is so rugged that it meets military shock and vibration standards. This analyzer requires no consumables and very little maintenance for low cost-of-ownership. Scientists using Picarro analyzers have reported the highest quality data, day in and day out, with fewer calibrations than other spectral absorption-based instruments.

Remote Operation: An Internet or dial-up modem connection can open up a world of possibilities for researchers. Users can connect remotely with the analyzer's internal Windows-based PC to control the instrument. The analyzer can automatically send data via email at regular intervals. The G2132-*i* synchronizes automatically with an atomic clock time service to timestamp the data accurately.

Performance Specifications		High Precision Mode	High Dynamic Range Mode	
Precision, δ^{13} C in CH ₄ (1- σ , 1 hr window)		< 0.8 % guaranteed precision at > 1.8 ppm 5 min. average < 0.5 % guaranteed precision at > 1.8 ppm, 15 minute average	< 0.4 % guaranteed precision at > 10 ppm	
Max Drift at STP δ ¹³ C in CH ₄ (over 24 hrs, peak-to-peak, 1 hr interval average)		< 1.5 % guaranteed drift at 10ppm < 2 % typical drift at 1.8ppm	< 1.5 ‰ guaranteed drift at 10 ppm	
Precision, CH ₄ Concentration (30 sec, 1-σ)		5 ppb + 0.05 % of reading (¹² C) 1 ppb + 0.05 % of reading (¹³ C)	50 ppb + 0.05 % of reading (¹² C) 10 ppb + 0.05 % of reading (¹³ C)	
Precision, CO ₂ Concentration (30 sec, 1-σ)		1 ppm + 0.25 % of reading (12C)		
Precision, H ₂ O Concentration (30 sec, 1-σ)		100 ppm		
CH₄ Dynamic Range		1.8-12 ppm guaranteed range 1.2-15 ppm operational range	10-1000 ppm guaranteed range 1.8-1500 ppm operational range	
CO₂ Dynamic Range		200 - 2000 ppm guaranteed range 0.01 - 0.4 % operational range		
H ₂ O Dynamic Range		0-2.4 % guaranteed range 0-5 % operational range		
Ambient Temperature Dependence		Guaranteed < ± 0.06 %/°C, typical < ± 0.025 %/°C		
Measurement Interval (includes periodic H ₂ O and CO ₂ measurement)		~ 2 secs	~ 2 secs	
se/Fall time (10-90% / 90-10%)		Typical ~30 secs		
Applications Considerations		Interference can occur for concentrations of H ₂ O and CO ₂ well above normal ambient levels, as well as other organics, ammonia, ethane, ethylene, or sulfur containing compounds. Users should verify with prepared lab samples. Please contact us to discuss the experimental conditions. Pressure drops in the instrument's gas path can draw external air when this system is used in recirculating applications.		
Analyzer Specifications				
Measurement Technique	CRI	CRDS		
Measurement Cell Temp. Control	+/- (+/- 0.005 °C		
Measurement Cell Pressure Control	+/- (+/- 0.0002 atm		
Shock and Vibration Testing	Meets shock and vibration MIL-STD 810F test standard and operates as specified afterward.			
Sample Temperature	-10 to 45 °C			
Sample Pressure	300 to 1000 Torr (40 to 133 kPa)			
Sample Flow Rate	< 50 sccm (typical ~25 sccm) at 760 Torr, no filtration required			
Sample Humidity	< 99 % RH non-condensing @ 40 °C, no drying required			
Ambient Temperature Range	10 to 35 °C (operating) -10 to 50 °C (storage)			
Ambient Humidity	< 99 % RH non-condensing			
Accessories	Pump (external), keyboard, mouse, LCD monitor (optional)			
Data Outputs	RS-232, Ethernet, USB, analog (optional) 0 – 10 V			
Fittings	1/4" Swagelok ®			
Dimensions	Analyzer: 17" w x 7" h x 17.55" d (43.18 x 17.78 x 44.57 cm) not incl. 0.5" feet External Pump: 5.6" w x 6.4" h x 11.9" d (14.3 x 16.3 x 30.3 cm)			
Installation	Ben	Benchtop or 19" rack mount chassis		
Weight	60.4	60.4 lbs (27.4 kg), includes external pump		

This product is not optimized for vehicular deployment where there is a requirement for pin-pointing precise methane source locations while driving. As a result, we do not support this product's use for natural gas leak detection or other real-time methane emissions applications while driving. The Picarro SurveyorTM system is the optimal product for such studies.

125 W (analyzer), 80 W (pump) at steady state

100 - 240 VAC, 47 - 63 Hz (auto-sensing), < 260 W start-up (total);

Power Requirements