

# Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>) Gas Concentration Analyzer

Continuous Monitoring for Sterilization  
& Bio-Decontamination Processes

# PICARRO



Avoid Oxidation of Sensitive Pharmaceuticals!

- Continuously measures parts-per-billion (ppb) levels of hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>)
- Built-in software ensures compliance with electronic records regulations
- Aeration cycle real-time monitoring
- Infrequent maintenance, no wet chemistry or consumables required

The **Picarro PI2124 gas concentration analyzer** ensures ultra-low residual hydrogen peroxide levels to help avoid oxidation and safeguard drug stability in isolated GMP pharmaceutical manufacturing applications. Production environment decontamination with vaporized hydrogen peroxide (VHP) is the industry standard. While VHP technology offers many performance and safety benefits, residual H<sub>2</sub>O<sub>2</sub> can easily oxidize pharmaceutical products. Biologic products are especially vulnerable, often showing significant oxidation even when exposed to low concentrations of H<sub>2</sub>O<sub>2</sub> (ppb range).

The Picarro PI2124 analyzer can continuously measure H<sub>2</sub>O<sub>2</sub> levels down to 3 ppb with precision better than 1 ppb. The analyzer is optimized for GMP applications. Installation and operation qualification is quick and easy,

as is validation with commercially available surrogate gases. The built-in software of the PI2124, designed to assist with 21 CFR part 11 compliance, offers strict regulations on electronic records, time stamps and audit trails. Combined with detailed validation reports, software and system documentation, the PI2124 offers comprehensive and tightly controlled data management software that complies with pharmaceutical and commercial manufacturing regulations.

Operating the PI2124 analyzer is simple and cost-effective—it requires no wet chemistry or consumables. Patented Picarro cavity ring-down spectroscopy (CRDS) technology has no moving parts and incorporates wavelength monitoring, providing long-term stability and reliability, as well as infrequent calibration and maintenance.

PI2124 Performance Specifications	Specifications**
<b>Precision (1σ, 10 sec)</b> <b>Precision (1σ, 300 sec)</b>	<3 ppb + 0.1% of reading <1 ppb + 0.1% of reading
<b>Lower Detection Limit (3σ, 300 sec)</b>	<3.0 ppb
<b>Zero Drift* (72 hours)</b> <b>Zero Drift (1 month)</b>	±2.5 ppb -5 ppb/+10 ppb
<b>Accuracy</b>	±5% of reading
<b>Measurement Range</b>	0–100 ppm
<b>Measurement Interval</b>	~10 sec
<b>Response Time</b> (Rise/Fall Time 10–90% / 90–10%)	<1 min

\* Picarro analyzers do not require a zero reference gas or zero cartridge to operate or meet specifications.

\*\* Specifications and an instrument-specific testing report (Certificate of Compliance) provided with every analyzer purchase.

## PI2124 System Specifications

Measurement Technique	Cavity Ring-Down Spectroscopy
Measurement Cell Temperature Control	±0.005°C
Measurement Cell Pressure Control	±0.0002 atm
Sample Temperature	-10 to 45°C
Sample Flow Rate	<1 slm at 760 Torr, no filtration required
Sample Pressure	300 to 1000 Torr (40 to 133 kPa)
Sample Humidity	<99% R.H. non-condensing @40°C, no drying required
Temperature	10 to 35°C (operating) -10 to 50°C (storage)
Ambient Humidity	<99% R.H. non-condensing
Other Gases Measured	H <sub>2</sub> O, CH <sub>4</sub>
Accessories	Pump (external, included), keyboard & mouse (included), LCD monitor (optional)
Communication Interfaces	RS-232, Ethernet, USB, analog 4-20mA, Modbus
Fittings	1/4" Swagelok® PFA Fittings
Dimensions	17" w x 8.38" h x 21" d (43.2 x 17.9 x 53.3 cm) including feet, not including small external pump module, 7.5" w x 4" h x 11" d (19 x 10.2 x 28 cm)
Weight	73 lbs (33.2 kg) including pump
Power Requirements	100–240 VAC, 50/60 Hz (auto-sensing), <260 W start up (total): 110 W (analyzer), 35 W (pump) at steady state
Certifications	CE Mark, CDRH
Country of Manufacture	USA
System Validation	Includes software wizard for analyzer validation using CH <sub>4</sub> as surrogate gas standard

## Software and User Accounts Regulations

Regulations Features	Not Signed In	Picarro Sign In	Operator	Technician	Administrator
<b>Config Menu/Settings</b> (requires additional sign in to access validation or user admin)		✓			✓
<b>View Data Viewer</b>	✓		✓	✓	✓
<b>Change Displayed Data</b>		✓	✓	✓	✓
<b>Quit Measuring</b>			✓	✓	✓
<b>Configure 4-20mA Output</b>		✓			
<b>H<sub>2</sub>O<sub>2</sub> Validation</b>				✓	✓
<b>Change displayed data</b> (full access)					✓
<b>User Administration</b>					✓

## PI2124 Validation Report

Picarro PI2124 Validation Report				
Operated by:	(username: admin)			
Start Time:	2025-11-04 09:45:10 (GMT+800)			
Signed by:	(username: admin)			
Signature Time:	2025-11-04 10:24:18 (GMT+800)			
<b>Summary</b>				
	Acceptance criteria	Result	Status	
Zero H <sub>2</sub> O <sub>2</sub>	>5ppb and < 10ppb	-1.556 ppb	Pass	
CH <sub>4</sub> slope	> 0.95 and < 1.05	0.999746	Pass	
CH <sub>4</sub> deviation	< 5%	0.343722%	Pass	
<b>Details</b>				
Step	Zero Air	Calibrant 1	Calibrant 2	Calibrant 3
Gas Source	cylinder1	cylinder2	cylinder3	cylinder4
Nominal CH <sub>4</sub> (ppm)	0.0±0.0%	7.056±0.9%	50.88±0.9%	100.4±0.8%
Observed average CH <sub>4</sub> (ppm)	0.011	7.080	50.822	100.404
CH <sub>4</sub> SD (ppm)	0.0634	0.0618	0.0643	0.0918
CH <sub>4</sub> deviation (%)	N/A	0.344	0.115	0.004
Observed H <sub>2</sub> O <sub>2</sub> (ppb)	-1.556	-0.461	-1.205	-2.105
Observed H <sub>2</sub> O (%)	0.001	0.002	0.000	0.000
CH <sub>4</sub> slope:	0.999746	CH <sub>4</sub> R <sup>2</sup> :	0.999999	
CH <sub>4</sub> intercept (ppm):	-0.005	Zero air CH <sub>4</sub> (ppm):	0.011	
H <sub>2</sub> O <sub>2</sub> equivalent (ppb):	-0.07	H <sub>2</sub> O <sub>2</sub> measured zero (ppb):	-1.556	

PICARRO

© 2026 PICARRO, INC.  
41-0099 Rev. B  
LIT CODE: 260203

3105 Patrick Henry Drive, Santa Clara, CA 95054 USA | 408-962-3900 | sales@picarro.com | picarro.com