PICARRO G2121-i Isotopic Carbon-Analyzer

Isotopic Carbon-Analyzer coupled to sample preparation systems

- · High-precision & high-accuracy with minimal drift
- Two-fold dynamic range (2000ppm-4000ppm)
- · Field and laboratory deployable with no consumables
- Installed and operational in 60 minutes
- · Rugged and insensitive to changes in ambient temperature
- · For use with CM-CRDS and the AutoMate Sample Preparation device only




Figure and data: Analysis of five acetanilide solid samples using CM-CRDS.

ΡΙΟΔ R R Ο

The World's Leading Instruments for Carbon and Water Cycle Measurements



Fig 1. Analysis of five acetanilide solid samples using Picarro Combustion Module (CM) coupled to the G2121-i CRDS. The standard deviation of the d13C values of the five sample replicates is outstanding.

Advantage Note: The Picarro G2121-*i* Isotopic CO₂ Analyzer is capable of measuring stable isotope ratios of carbon in CO₂ and total CO₂ concentration when coupled to a sample preparation front-end such as an Elemental Analyzer (EA), a Dissolved Inorganic Carbon acidification module or with Picarro's Combustion Module (CM) designed by Costech. The analyzer is based on Picarro's unique Cavity Ring Down Spectroscopy (CRDS), a time-based measurement utilizing a near-infrared laser to measure a spectral signature of the molecule. Gas is circulated in an optical measurement cavity with an effective path length of up to 20 kilometers. A patented, high-precision wavelength monitor makes certain that only the spectral feature of interest is being monitored, greatly reducing the analyzer's sensitivity to interfering gas species, and enabling the high-precision ₁₃C isotope measurements even if there are other gases present. As a result, the analyzer maintains high linearity, high precision, and high accuracy over a two-fold dynamic range (2000ppm-4000ppm) with minimal calibration required. Precise temperature and pressure control systems designed into the Picarro G2121-*i* ensure accurate measurements over long periods of time with minimal use of calibration gases or standard reference material. The analyzer itself is exceptionally rugged and insensitive to changing ambient temperature, essentially drift and maintenance free, and requires no consumables, thereby offering significant ease-of-use and cost of ownership benefits.

Easily transportable from site to site, the analyzer can be set up and acquiring data within 60 minutes. The gas concentration is displayed in real-time with no post-processing required, and is continuously archived to the analyzer's internal hard drive. Designed to operate both in laboratories and in field settings, the analyzer can be remotely accessed for data monitoring with the analyzer's internal Windows-based PC and control it through a standard

Remote Desktop connection or with similar remote login software. The analyzer can also use the internet connection to automatically synchronize with an atomic clock time service.

ΡΙСΔ R R O

The World's Leading Instruments for Carbon and Water Cycle Measurements

Note: If customers wish to use the G2121-*i* for atmospheric applications, they will need to purchase the Picarro Ambient Air Certification (P/N A0502) or they may

expect a significant reduction in precision and accuracy at ambient air conditions. Also, it is important to remember that in the case of dual use (ambient and high concentration) two separate calibrations should be performed; one for the high concentration CO₂ level with N₂ as a background gas and one at ambient air conditions. Picarro provides an isotopic calibration software utility that streamlines this process.

About Picarro CRDS

Picarro's unique Cavity Ring-Down Spectroscopy (CRDS) system is a time-based measurement system that uses a near-infrared laser to quantify spectral features of molecules in a sample gas passed through an optical measurement cavity. An effective path length of up to 20 kilometers inside the cavity results in exceptional instrument precision and sensitivity. The optical cavity is thermally controlled to ensure minimal drift, even in the harshest environments. Precise temperature and pressure control systems designed into the analyzer ensure accurate measurements over long periods of time as well as allow extremely infrequent calibration. Picarro CRDS systems all include a patented, high-precision wavelength monitor that makes certain the analyzer only monitors spectral features of interest. For researchers, the CM-CRDS delivers a best-in-class combination of flexibility, speed, high-precision, and ease of use that sets a new standard for bulk stable isotope analysis.

Targeted Performance Specification	
Gas Species	Precision (5min pulse average)
δ ¹³ C in CO ²	<0.4 ‰ (0.2-0.3 ‰ typical) sample-to-sample
12CO2 Concentration Range	2000ppm – 4000ppm

System Specifications	
Measurement Technique	CRDS
Range	CO ₂ : 2000-4000 ppm
Measurement Interval	~10 seconds
Sample Temperature	-10 to 45 °C
Sample Flow Rate	< 0.1 slm at 760 Torr, no filtration required
Sample Pressure	300 to 1000 Torr (40 to 133 kPa)
Sample Humidity	< 99% R.H. non-condensing @ 40 °C, no drying required
Temperature	10 to 35 °C (operating) -10 to 50 °C (storage)
Humidity (ambient)	< 99% R.H. non-condensing
Other Gases Measured	H ₂ O
Accessories	Pump (external), keyboard, mouse, LCD monitor (optional)
Outputs	RS-232, Ethernet, USB, analog (optional) 4-20 mA / -10 – 10V
Fittings	1/4" Swagelok ®
Dimensions	Analyzer: 17" w x 7" h x 17.55" d (43.18 x 17.78 x 44.57 cm) not incl. 0.5" feet External Pump: 5.6" w x 6.4" h x 11.9" d (14.3 x 16.3 x 30.3 cm)
Installation	Benchtop
Weight	60.4 lbs (27.4kg), includes external pump
Power Requirements	90-120 VAC, 50/60 Hz, 220 VAC, 50 Hz, < 300 W