

# PICARRO

# S3110

## Quick Start Guide

Yearly Maintenance Kit for G4xxx Analyzers

### Maintenance Kit Box Contents:

#### Particulate filter, Teflon (PN S1021)

User-replaceable particulate filter for removing particulates from sample stream. Clogged filter will reduce sample flow.



*This maintenance kit is compatible with the following analyzers:*

*G4301, G4302*

# PICARRO

# S3110

## Quick Start Guide

Yearly Maintenance Kit for G4xxx Analyzers

### Fan and Filter Maintenance and Replacement Instructions:

For detailed instructions on how to replace filters and fans, please refer to the user manuals and how-to videos listed below.

#### **Picarro Video Tutorial:**

##### **Replacing Your Inlet Particulate Filter**

Inlet particulate filters on Picarro instruments should be changed every year. This tutorial provides step-by-step instructions for replacing and checking your inlet filter for leaks.

[vimeo.com/375518688](https://vimeo.com/375518688)

#### **Detailed Maintenance Guide:**

##### **Maintenance Guide: Inlet Particulate Filter Replacement for Picarro Analyzers**

(Document Number 40-0086)

[bit.ly/40-0086](https://bit.ly/40-0086)

41-0067 Rev B

### Need technical support?

Go online: Visit [picarro.com/support](https://picarro.com/support)

Email us: [support@picarro.com](mailto:support@picarro.com)

Call us: +1 408-962-3991 (USA center) +31 858-881-650 (European center)